



Lohmann & Rauscher

Compression Stockings and Garments



Simple. Effective. Swiss quality.



Return to
contents page

VENOSAN® Compression Stockings help to make life easier in many situations. The pressure VENOSAN Compression Stockings exert promotes healthy veins. They can help to prevent venous disease as well as being an important part of the treatment for patients being treated for acute and chronic venous disorders.

The diverse VENOSAN range offers the right product for a wide variety of needs. All VENOSAN Compression Stockings have the following characteristics in common:

Simple.

The VENOSAN range provides a simple guide to help you determine the correct size, compression class and style of compression stocking.

The wear and care instructions included with VENOSAN Compression Stockings are easy to follow which patients appreciate.

Effective.

VENOSAN Compression Stockings are RAL-certified, meaning that they provide effective therapeutic compression. The effective design of VENOSAN Compression Stockings means that for patients, many will experience enhanced treatment success and improved quality of life.

Swiss quality.

Products from Switzerland are known for excellent workmanship and top quality.

Manufactured using high quality materials, such as the original Lycra® and elastane, VENOSAN Compression Stockings offer wearing comfort and durability.



Product lines

Medical compression stockings and pantyhose.



Return to contents page

| Line | Colours | CCL | | | | | | | | |
|---|----------------|-------|------------------------------|--|------------------------------------|--|--|--|-----------------|-------------------------------|
| | | | AD Below knee stocking | ADH Below knee stocking with dotted top band | AG Thigh- length stocking | AGH Thigh-length stocking with dotted top band | AGH Thigh-length stocking with lace top band | AGG Thigh- length stocking with belt | AT Pantyhose | ATU Maternity pantyhose |
| VENOSAN® 3000 NATURAL LOOK | Beige | CCL 2 | ● | ● | | ● | | | | |
| | | CCL 3 | ● | | | | | | | |
| VENOSAN® 4000 CLIMATE EFFECT | Mexico | CCL 1 | ● | ● | ○ | ● | ● | ○ | ● | ● |
| | Marokko | CCL 2 | ● | ● | ○ | ● | ● | ● | ● | ● |
| | Black | CCL 2 | ● | ● | ○ | ● | ● | ● | ● | ● |
| VENOSAN® 5000 THE MAGIC OF SILVER | Mexico | CCL 1 | ● | ● | ○ | ● | ● | ● | ● | ○ |
| | Marokko | CCL 2 | ● | ● | ○ | ● | ● | ● | ● | ○ |
| | Black | CCL 3 | ● | ○ | ○ | ○ | ○ | ○ | ○ | ○ |
| VENOSAN® 6000 CLASSIC | Beige Black | CCL 1 | ● | | | ● | | | ● | |
| | | CCL 2 | ● | | | ● | | ● | ● | |
| | | CCL 3 | ● | | | ● | | | ● | |
| VENOSAN® 7000 NATURAL TOUCH | Beige Black | CCL 1 | ● | | | | | | | |
| | | CCL 2 | ● | ● | ○ | ● | ○ | ● | ○ | |
| | | CCL 3 | ● | ● | ○ | ● | ○ | ○ | ○ | |
| VENOSAN® 8000 ULCERFIT SUCCESSFUL THERAPY | Beige | CCL 2 | ● | | | | | | | |
| | | CCL 3 | ● | | | | | | | |
| VENOSAN® SOFT VENOSAN® PRO VENOSAN® STRONG FLAT KNIT THERAPY | Beige Black | CCL 1 | ○ | ○ | ○ | ○ | | ○ | ○ | |
| | | CCL 2 | ○ | ○ | ○ | ○ | | ○ | ○ | |
| | | CCL 3 | ○ | ○ | ○ | ○ | | ○ | ○ | |
| | | CCL 4 | ○ | ○ | ○ | ○ | | ○ | ○ | |

● available in standard sizes ○ available in custom sizes only ● available in standard and custom sizes

Pressure at ankle **CCL 1** 18–21mmHg **CCL 2** 23–32mmHg **CCL 3** 34–46mmHg **CCL 4** >49mmHg

Product overview

Medical compression stockings and pantyhose.




Return to
contents page

VENOSAN® 3000 COTTON NATURAL LOOK

- Durable
- 100% cotton on the skin
- Absorbs moisture
- Easy to don

Indications

- CVI (chronic venous insufficiency) 2–3
- Lipoedema stage 2–3
- After vein surgery
- PTS (post-thrombotic syndrome)
- Mild lymphoedema



Composition


- CCL 2 49% cotton
27% nylon polyamide
24% elastane
- CCL 3 40% cotton
27% nylon polyamide
33% elastane

VENOSAN® 4000 CLIMATE EFFECT

- Climate effect
- Semi-transparent
- Comfortable
- Highly elastic

Indications

- CVI (chronic venous insufficiency) 1
- Varicosis during pregnancy
- After vein surgery



Composition


- 71% TACTEL® polyamide
29% elastane

VENOSAN® 5000 THE MAGIC OF SILVER

- Odour-inhibiting
- Anti-static
- Temperature-regulation
- Anti-bacterial

Indications

- CVI (chronic venous insufficiency) 1–2
- After vein surgery
- Mild skin conditions
- Lipoedema stage 1



Composition

- 40% TACTEL® polyamide
30% elastane
20% cotton
10% polyamide/silver (X-Static®)

Product overview

Medical compression stockings and pantyhose.



Return to contents page

VENOSAN® 6000 CLASSIC

- Maximum sheerness
- Breathable
- Ease of care

Indications

- CVI (chronic venous insufficiency) 1–2
- Varicose veins
- Varicosis during pregnancy
- After vein surgery



Composition

70% nylon polyamide
30% elastane

VENOSAN® 7000 NATURAL TOUCH

- High stiffness (resistance)
- Skin-friendly
- Temperature regulation

Indications

- CVI (chronic venous insufficiency) 2–3
- Lipoedema stage 2–3
- After vein surgery
- After thrombosis
- PTS (post-thrombotic syndrome)
- Mild lymphoedema
- Marked orthostatic stress



Composition


62% TACTEL® polyamide
26% elastane
8% cotton
4% Lyocell/SeaCell®

VENOSAN® 8000 ULCERFIT SUCCESSFUL THERAPY

- Two-layer stocking system
- Fast and easy to use
- High wearing comfort
- Liner can be worn at night

Indications

- Treatment of venous leg ulcers
- Helps to prevent recurrence of venous leg ulcers



Composition

Overstocking:
62% TACTEL® Polyamide
26% Elastane
8% Cotton
4% Lyocell/SeaCell®

Liner:
90% Polyamide
10% Elastane

VENOSAN® SOFT VENOSAN® PRO VENOSAN® STRONG FLAT KNIT THERAPY

- Flat-knit custom made garments
- Smooth and durable
- Low to high stiffness
- Manufactured to European standards

Indications

- Treatment and management of Chronic Oedema and Lymphoedema



Composition

Polyamide & Elastane
VENOSAN PRO (CCL II)
& VENOSAN STRONG (CCL 4S)
also contains cotton
VENOSAN PRO also contains Latex



Return to contents page

VENOSAN® 3000 COTTON

NATURAL LOOK

100% cotton - natural and soft.



■ Natural

Cotton allows a higher degree of wearer comfort for all those who prefer a natural fibre against the skin, or who are allergic to synthetic fibres. 100% cotton against the skin.

■ Durable

The construction ensures the compression stocking is durable and hard wearing making it ideal for long-term compression therapy.


■ Easy to don

The stretchable knitted cotton fabric allows easy donning of the compression stocking.

[Click for VENOSAN 3000 ordering and sizing details](#)

Compression Stockings and Garments

Beige



CCL 2

CCL 3

Garment Styles



Indications

- CVI 2–3 (chronic venous insufficiency)
- Lipoedema stage 2–3
- Post vein surgery
- PTS (post-thrombotic syndrome)
- Mild lymphoedema

Composition

- CCL 2: 49% cotton, 27% nylon polyamide, 24% elastane
- CCL 3: 40% cotton, 27% nylon polyamide, 33% elastane



Return to contents page

VENOSAN® 4000

CLIMATE EFFECT

Microfibre - fine and stylish.



■ Climate effect

The air permeable yarns and knitting construction of this stocking allow for moisture to be wicked away from the skin and through the stocking.

■ Fashionably elegant

Wearable in all seasons. It has a contemporary stylish look, with the sheer transparency.



■ Comfortable

The velvety touch of the high tech TACTEL® microfibrils offers outstanding comfort. The microfibrils are many times thinner than a human hair.

■ Made to Measure available

[Click for VENOSAN 4000 ordering and sizing details](#)

Compression Stockings and Garments



Garment Styles



Indications

- CVI 1 (chronic venous insufficiency)
- Pregnancy-induced varicose veins
- Milder symptoms
- Post vein surgery

Composition

- 71% TACTEL® polyamide
- 29% elastane



Return to contents page

VENOSAN® 5000

THE MAGIC OF SILVER

For ideal temperature regulation.



■ Durable

The innovative yarn combination of silver and cotton make this stocking suitable for active men and women.

■ Regulates temperature



Silver is well known for its excellent heat conductive properties. In combination with the microfibre TACTEL®, silver produces rapid dissipation of skin moisture, helping to ensure optimal wearing comfort.

■ Silver-coated X-Static® fibres

Contains silver-coated X-Static fibre. No reduction of the properties of silver even after 150 washing cycles.

■ Made to Measure available

[Click for VENOSAN 5000 ordering and sizing details](#)



Garment Styles



Indications

- CVI 1–2 (chronic venous insufficiency)
- Post vein surgery
- Mild skin conditions
- Lipoedema stage 1

Composition

- 40% TACTEL® polyamide
- 30% elastane
- 20% cotton
- 10% polyamide / silver X-Static®



Return to contents page

VENOSAN® 6000

CLASSIC


Reliable and cost-effective.



- **Durable knit construction**
Basic yarns make the stocking cost-effective.
- **Anatomically knitted**
For maximum comfort, does not slip.
- **Ease of care**
Can be machine washed at 40°C on a delicates programme. Dries naturally overnight.

[Click for VENOSAN 6000 ordering and sizing details](#)

Compression Stockings and Garments



Garment Styles



Indications

- CVI 1-2 (chronic venous insufficiency)
- Varicose veins
- Varicosis during pregnancy
- After vein surgery

Composition

- 70% nylon polyamide
- 30% elastane



Return to contents page

VENOSAN® 7000

NATURAL TOUCH

High stiffness with optimal comfort.



■ High stiffness

Due to the enhanced stiffness the stockings provide a massage effect and are characterised as shaping, opaque compression stockings with high stiffness.

■ Skin-friendly

The composition ensures both high resistance and optimal wearing comfort.


■ Strong natural fibres

Due to the durability of the material these stockings are suitable for people who demand a lot from their stockings at work, in sports and in their leisure time.

■ Made to Measure available

[Click for VENOSAN 7000 ordering and sizing details](#)

Compression Stockings and Garments



Garment Styles



Indications

- CVI 2–3 (chronic venous insufficiency)
- Lipoedema stage 2–3
- Following acute thrombosis (following dressing)
- Post vein surgery
- PTS (post-thrombotic syndrome)
- Mild lymphoedema
- Marked orthostatic stress

Composition

- 62% TACTEL® polyamide
- 26% elastane
- 8% cotton
- 4% Lyocell/SeaCell® (seaweed fibre)



Return to contents page

VENOSAN® 8000 ULCERFIT

SUCCESSFUL THERAPY

Compression Stocking System.




- **Two layer system**
Offers optimum therapy thanks to the ideal combination of liner and overstocking.
- **Time saving treatment**
Fast and easy to use, patient can apply the system without assistance.
- **Better patient compliance**
High wearer comfort. The flat and thin surface of the liner ensures easy donning and doffing of the overstocking.

[Click for VENOSAN 8000 ordering and sizing details](#)

Compression Stockings and Garments

Beige



CCL 2 CCL 3

Version



Indications

- Treatment of venous leg ulcers
- Helps to prevent recurrence of venous leg ulcers

Composition

Overstocking:

- 62% TACTEL® Polyamide
- 26% Elastane
- 8% Cotton
- 4% Lyocell/SeaCell®

Liner:

- 90% Polyamide
- 10% Elastane



Return to contents page

VENOSAN® 4000 Armsleeves

CLIMATE EFFECT

Microfibre - fine and stylish.



■ Climate effect

The air permeable yarns and knitting construction of this stocking allow for moisture to be wicked away from the skin and through the stocking.

■ Fashionably elegant

Wearable in all seasons. It has a contemporary stylish look, with the sheer transparency.

■ Comfortable


The velvety touch of the high tech TACTEL® microfibrils offers outstanding comfort. The microfibrils are many times thinner than a human hair.

■ Made to Measure available

[Click for VENOSAN 4000 Armsleeves ordering and sizing details](#)


Compression Stockings and Garments

Beige



CCL 2

Garment Styles



Indications

- Mild lymphoedema

Composition

- 71% TACTEL® polyamide
- 29% elastane



Return to contents page

VENOSAN® 7000 Armsleeves

NATURAL TOUCH

High stiffness with optimal comfort.




- **High stiffness**
Due to the enhanced stiffness the stockings provide a massage effect and are characterised as shaping, opaque compression stockings with high stiffness.
- **Skin-friendly**
The composition ensures both high resistance and optimal wearing comfort.
- **Strong natural fibres**
Due to the durability of the material these stockings are suitable for people who demand a lot from their stockings at work, in sports and in their leisure time.
- **Made to Measure available**


[Click for VENOSAN 7000 Armsleeves ordering and sizing details](#)

Compression Stockings and Garments

Beige



Garment Styles



Indications

- Mild lymphoedema

Composition

- 62% TACTEL® polyamide
- 26% elastane
- 8% cotton
- 4 % Lyocell/SeaCell® (seaweed fibre)