


NOSTRIL RETAINER [non-sterile]




■ Features

1. This product prevents the deformation of the nasal foramen after the surgery for cleft lip or palate, rhinoplasty, or augmentation rhinoplasty without imposing a burden on the patient.
2. 13 sizes each for the standard and long types are available so that the appropriate size for each patient can be selected.
3. Made of soft, less irritative silicone.

● Nostril Retainer



● Mold Guide



■ Precautions for use

- Do not damage the product while picking with forceps.
The product may be torn from the damaged part.
- Once it is damaged, please use new product.

● Sterilization Method: Autoclaving

Standard Type

| Product No. | Size W × L × t (mm) | Packing |
|-------------|---------------------|---------|
| #8101 | No.1 16 × 8 × 7 | 1pc/box |
| #8102 | No.2 17 × 9 × 8 | 1pc/box |
| #8103 | No.3 18 × 10 × 8 | 1pc/box |
| #8104 | No.4 19 × 10 × 9 | 1pc/box |
| #8105 | No.5 20 × 10 × 10 | 1pc/box |
| #8106 | No.6 21 × 12 × 11 | 1pc/box |
| #8107 | No.7 22 × 12 × 11 | 1pc/box |
| #8108 | No.8 24 × 13 × 13 | 1pc/box |
| #8109 | No.9 25 × 14 × 14 | 1pc/box |
| #8110 | No.10 26 × 15 × 15 | 1pc/box |
| #8111 | No.11 27 × 16 × 16 | 1pc/box |
| #8112 | No.12 29 × 17 × 17 | 1pc/box |
| #8113 | No.13 30 × 17 × 18 | 1pc/box |

Long Type

| Product No. | Size W × L × t (mm) | Packing |
|-------------|---------------------|---------|
| #8131 | No.1 16 × 14 × 7 | 1pc/box |
| #8132 | No.2 17 × 15 × 8 | 1pc/box |
| #8133 | No.3 18 × 16 × 8 | 1pc/box |
| #8134 | No.4 19 × 16 × 9 | 1pc/box |
| #8135 | No.5 20 × 16 × 10 | 1pc/box |
| #8136 | No.6 21 × 18 × 11 | 1pc/box |
| #8137 | No.7 22 × 18 × 11 | 1pc/box |
| #8138 | No.8 24 × 19 × 13 | 1pc/box |
| #8139 | No.9 25 × 20 × 14 | 1pc/box |
| #8140 | No.10 26 × 21 × 15 | 1pc/box |
| #8141 | No.11 27 × 22 × 16 | 1pc/box |
| #8142 | No.12 29 × 23 × 17 | 1pc/box |
| #8143 | No.13 30 × 23 × 18 | 1pc/box |

Mold Guide

| Product No. | Product Name | Packing |
|-------------|------------------------------|----------|
| #8100 | Mold Guide for Standard Type | 1set/box |
| #8150 | Mold Guide for Long Type | 1set/box |

* Mold Guide for proper size selection.