

Siltape[®]

Soft silicone perforated tape

Preventing facial skin damage beneath PPE with Siltape

During the COVID-19 outbreak, healthcare professionals are seeing early signs of skin damage with the increased use of facial masks. Accurate fit testing and the use of barrier creams can help. If additional protection is needed under a mask, Siltape can be used.*

*When using fit tested masks, ensure the fit test is done after the application of Siltape and that the seal isn't compromised.

What is Siltape?

Siltape is a soft silicone tape on a handy roll with perforations for ease of use; simply tear off what you need and apply to the affected area.

Advancis
Medical

Where to use Siltape:

- Bridge of nose
- Cheeks
- Tops of ears
- Chin
- Forehead



Ordering Information

Size	Stock code	Pack size
2cm x 3m	CR3938	1
4cm x 1.5m	CR3939	1



for further information please contact us:

P 0800 226 3673 W bamford.co.nz E enquiries@bamford.co.nz

12 Victoria St, Lower Hutt, 5010, New Zealand