



Transparent Dressings



OVERVIEW

Transparent film wound dressing, allowing visual inspection of the wound site whilst promoting a “moist wound healing environment”. Newly forming tissue is protected whilst a wound is protected from contamination by dirt or bacteria.

- Gentle self-adhesive dressing
- Can be used in the shower
- High moisture vapour transfer rate (MVTR)

CODE	DESCRIPTION	STERILE	SELLING UNIT	OUTER QTY
36590006	Transparent Dressing 4cm x 5cm	Sterile	Box/100	20
36590007	Transparent Dressing 6cm x 7cm	Sterile	Box/100	20
36590008	Transparent Dressing 10cm x 12cm	Sterile	Box/50	16
36590010	Transparent Dressing 10cm x 25cm	Sterile	Box/30	10
36590051	Transparent Dressing 15cm x 20cm	Sterile	Box/25	10
36590052	Transparent Dressing 10cm x 14cm	Sterile	Box/50	16