



Non-Woven Island Dressings



OVERVIEW

Non-Woven Island dressings are made with a breathable non woven backing which is absorbent with a non adherent pad. Providing patient protection and comfort offering Secure and reliable fixation which Conforms to skin contours.

Ideal for use on

- Post operative wounds
- In A&E for management of cuts and lacerations
- For wounds with moderate exudate

Individually wrapped, sterile

CODE	DESCRIPTION	STERILE	SELLING UNIT	OUTER QTY
36590000	Non-Woven Island Dressing 8cm x 6cm	Sterile	Box/50	TBA
36590001	Non-Woven Island Dressing 5cm x 7.2cm	Sterile	Box/50	60
36590002	Non-Woven Island Dressing 10cm x 10cm	Sterile	Box/50	24
36590003	Non-Woven Island Dressing 10cm x 15cm	Sterile	Box/50	12
36590004	Non-Woven Island Dressing 10cm x 20cm	Sterile	Box/25	12
36590005	Non-Woven Island Dressing 10cm x 25cm	Sterile	Box/25	12
36590016	Non-Woven Island Dressing 8cm x 15cm	Sterile	Box/50	12
36590017	Non-Woven Island Dressing 8cm x 10cm	Sterile	Box/50	12
36590039	Non-Woven Island Dressing 6cm x 8cm	Sterile	Box/60	60