



Transparent Island Dressings



OVERVIEW

Transparent Island Dressings are sterile, self adhesive, transparent film wound dressings incorporating an absorbent, non adherent pad.

Made from a latex free semi permeable polyurethane film membrane coated with an hypoallergenic acrylic adhesive, the dressings are breathable, allowing good oxygen and moisture vapour exchange.

It is showerproof and impermeable to liquids, bacteria, and viruses.

Transparent Island Dressings are primarily indicated for use on acute wounds with low levels of exudate.

- Gentle self-adhesive dressing
- Clear, showerproof film
- Absorbent pad
- High moisture vapour transfer rate (MVTR)

CODE	DESCRIPTION	STERILE	SELLING UNIT	OUTER QTY
36590019	Transparent Island Dressing 5cm x 7.2cm	Sterile	Box/50	90
36590021	Transparent Island Dressing 12cm x 10cm	Sterile	Box/50	24
36590022	Transparent Island Dressing 8.5cm x 15.5cm	Sterile	Box/50	12
36590023	Transparent Island Dressing 20cm x 10cm	Sterile	Box/50	12
36590024	Transparent Island Dressing 10cm x 25cm	Sterile	Box/50	12
36590025	Transparent Island Dressing 30cm x 10cm	Sterile	Box/50	12
36590026	Transparent Island Dressing 8.5cm x 9.5cm	Sterile	Box/50	12
36590027	Transparent Island Dressing 35cm x 10cm	Sterile	Box/50	50