

**VENOSAN®**

# LEGLINE 20

Support Stockings 140 den

**ELEGANT LUXURY**



**VENOSAN®**

# LEGLINE 20

20 mm Hg  
140 den

Support Stockings 140 den

HEEL AND TOE REINFORCEMENT



## Stützstrümpfe | Support stockings | Bas de support | Calze riposanti

**AD**  
Wadenstrümpfe  
Below knee  
Bas au mollet  
Gambaletti

**AGH**  
Schenkelstrümpfe  
Thigh length  
Bas à cuisse  
Calza alla coscia

**AT\***  
Strumpfhose  
Pantyhose/Tights  
Collants  
Collant

**ATU Materna\***  
Schwangerschaftsstrumpfhose  
Maternity pantyhose  
Collants de grossesse  
Collant gravidanzaa

\* Komfort-Hosenteil  
Comfort panty  
Culotte confort  
Corpino comfort

Stufenlos abnehmende Kompression | Graduated compression |  
Compression dégressive | Compressione graduata

Umfang | Circumference | Circonférence | Circonferenza (cm)

AD	AGH, AT, ATU	
	60 % 12 mm Hg	cF
80 % 16 mm Hg	80 % 16 mm Hg	cD
		cC
100 % 20 mm Hg	100 % 20 mm Hg	cB

XS	S	M	L	XL
42-44	45-47	48-50	51-53	54-60
29-31	32-34	35-37	38-40	41-43
31-33	34-36	37-39	40-42	43-45
20-21	22-23	24-25	26-27	28-29

**SWISS COMFORT**

### Effective Compression

Verhindert geschwollene Beine  
Prevents swollen legs  
Évite des jambes enflées  
Impedisce le gambe gonfie

### Everyday Elegance

Transparente Optik  
Transparent look  
Look transparent  
Look trasparente

### Optimal Comfort

Für höchsten Tragekomfort, rutscht nicht  
For maximum comfort, does not slip  
Pour un confort optimal, ne glisse pas  
Non stringe, non scivola

Spitzenhaftband | Lace top band | Adhésive ornée | Bordo autoreggente pizzo




Aussen  
Outside  
Face externe  
Lato esterno



Innen  
Inside  
Côte intérieur  
Lato interno

Grösse | Height | Grandeur | Altezza



Artikel | Article | Article | Articolo

Type	*Colours	XS	S	M	L	XL
AD	Nude	VL 303 NU	VL 304 NU	VL 305 NU	VL 306 NU	-
	Sahara	VL 303 SA	VL 304 SA	VL 305 SA	VL 306 SA	-
	Black	VL 303 BL	VL 304 BL	VL 305 BL	VL 306 BL	-
AGH	Nude	VL 323 NU	VL 324 NU	VL 325 NU	VL 326 NU	-
	Sahara	VL 323 SA	VL 324 SA	VL 325 SA	VL 326 SA	-
	Black	VL 323 BL	VL 324 BL	VL 325 BL	VL 326 BL	-
AT	Nude	VL 333 NU	VL 334 NU	VL 335 NU	VL 336 NU	VL 337 NU
	Sahara	VL 333 SA	VL 334 SA	VL 335 SA	VL 336 SA	VL 337 SA
	Black	VL 333 BL	VL 334 BL	VL 335 BL	VL 336 BL	VL 337 BL
ATU	Sahara	VL 343 SA	VL 344 SA	VL 345 SA	VL 346 SA	VL347 SA
	Black	VL 343 BL	VL 344 BL	VL 345 BL	VL 346 BL	VL 347 BL

\*Colours



Nude



Sahara



Black



Materialzusammensetzung | Composition | Composition | Composizione

78 % Polyamide | 22 % LYCRA® Elastane





**VENOSAN®**  
Compression Stockings

[www.venosan.com](http://www.venosan.com)

VENOSAN® is a registered  
trademark of SWISSLASTIC AG  
ST. GALLEN, Switzerland

an  Company



SWISSLASTIC AG ST. GALLEN  
Rorschacher Strasse 304  
CH-9016 St. Gallen  
Switzerland