


VENOSAN® Support Stockings Product guide





Product overview Medical support stockings.



VENOSAN[®] MICROFIBERLINE STYLISH COMFORT, FASHIONABLY SOFT

- Soft and comfortable due to microfibre material
- Stylish support
- Non-restricting top band
- Reinforced heel and toe
- Available in men's and ladies' variants

Indications

- Tired, heavy, swollen and aching legs
- Prevention of Deep Vein Thrombosis (DVT)





AD

Below knee

stocking



VENOSAN[®] SILVERLINE EVERYDAY COMFORT

- Enhanced comfort and fit with knitted top band
- Elimination of bacteria and fungus associated with foot odour
- Anti-static effect
- Soft, skin friendly material with ribbed texture
- Wicks away moisture to keep the skin dry and comfortable
- Available in men's and ladies' variants

Indications

- Tired, heavy, swollen and aching legs
- Prevention of Deep Vein Thrombosis (DVT)



15-20mmHg

Composition

Men: 86% Nylon 9% Spandex 5% X-STATIC® Silver plated Nylon Ladies: 80% Nylon 11% Spandex 9% X-STATIC® Silver plated Nylon

Men

Black



Ladies

Black

Beige

VENOSAN[®] ATHLETIC PERFORMANCE SUPPORT

- Unisex compression socks
- Reciprocated heel and toe for improved fit & comfort
- Soft, skin-friendly material with ribbed texture
- Ventilated foot and heel for breathability
- Cushioned achilles tendon area for better protection and comfort

Indications

- Tired, heavy, swollen and aching legs
- Prevention of Deep Vein Thrombosis (DVT)



- Unisex compression socks
- Made with COOLMAX[®] fibres that provide superior wicking to keep feet and legs dry
- Cushioned comfort sole with extra padding support
- Wide toe box for roomy comfort
- Easy to put on

Indications

- Tired, heavy, swollen and aching legs
- Prevention of Deep Vein Thrombosis (DVT)



White

20-30mmHg

20-30mmHg

Composition 75% Nylon 12% Polyester

13% Spandex

White

Black

Pink

AD

Below knee

stocking

AD

Below knee

stocking

15-20mmHg



Blue

Purple


Anti-embolism stockings



- Unisex compression socks
- Anatomical design providing a good fit
- Colour coding for simple and reliable
- identification of size
- Top inspection window for safety

Indications

- To prevent thrombosis and embolism before, during and after surgeries
- General immobility, confinement to bed
- Before, during and after childbirth





18mmHg

Black

91% hylon polyamide 9% LYCRA® elastane AG: 89% hylon polyamide 11% LYCRA® elastane

White



Accessories for compression stockings.



VENO-GLIDER

Foot Slip Donning Aid

Article no.	Description	Unit
995 345	VENO-GLIDER	1 pc.



VENOSAN[®] Dot Gloves

Washable gloves for donning and doffing compression stockings, available in two sizes.

Material composition

Knitted fabric: 100% nylon polyamide, Dots: PVC

Article no.	Description	Unit
VC 003	VENOSAN Dot Gloves Size S/M	1 pair
VC 004	VENOSAN Dot Gloves Size L/XL	1 pair



COLTEX

Special adhesive for securing medical compression stockings to the skin.

Composition

PVP/VA copolymer, glycerin, 2-propanol

Article no.	Description	Unit
KC 67000	COLTEX Adhesive 60ml	1 pc.



VENOSAN® Eureka On

Donning aid for compression stockings. Can be used with open and closed toe.

Article no.	Description	Unit
155380	VENOSAN Eureka On	1 pc.

Lohmann & Rauscher Pty Ltd · 321 Ferntree Gully Road · Mount Waverley, Victoria 3149, Australia · Customer Service: 1300 572 869 Lohmann & Rauscher Limited · 21 Queen Street · Auckland 1010 · New Zealand · Customer Service: 0800 572 869 VENOSAN[®] is a registered trademark of SWISSLASTIC AG St. Gallen, Switzerland – an 認 Company COOLMAX[®] and LYCRA[®] are trademarks of The LYCRA Company