



WM Bamford & Co Ltd  
12 Victoria Street  
Private Bag 31346  
Lower Hutt 5040  
New Zealand  
tel +64 4 576 2100  
fax +64 4 569 6489  
enquiries@bamford.co.nz  
www.bamford.co.nz



## Compact Surgical LED Loupe Light with S-Frame eye protection (magnification Loupe lenses optional)

### SEEING IS BELIEVING

**Advanced LED Technology.** White, homogenous light, ideal for all surgical applications.

**Coaxial Illumination.** Shadow-free illumination of cavities, ideal. Ideal co-axial viewing with all HR/HRP loupes.

**Homogenous illumination of the whole field of view.** The light spot is ideal for HR / HRP loupe systems. 81mm illumination field at 420mm, 95mm at 500mm working distance.

**Compact, lightweight (24g).** Ensures perfect comfort.

**Aluminium housing.** Robust, long lasting.

**Waterproof design.** Easy to clean.

**Fully-portable with mPack LL.** Autonomy at least 17 hours.

**Individual adjustment of the angle of view.** The right angle for every situation.

**Registered Design.**

**40.000 Lux / 25 cm distance**

**No light box or Fibre Optic cables required**

When performance counts Bamford delivers