



Lohmann & Rauscher

ReadyWrap™

Clinically effective¹ and easy-to-use^{2,3}
adjustable compression wrap



ReadyWrap

Adjustable compression garments for patients with venous and lymphatic diseases.

ReadyWrap

Clinically-effective and easy-to-use.

ReadyWrap is an adjustable compression wrap that provides patients and their caregivers an easy-to-use solution to support wound healing and oedema management.¹ ReadyWrap provides clinically effective compression therapy while supporting patient compliance with its unique alternating strapping system.



“To manage the swelling reduction, I placed him in bi-lateral ReadyWrap garments. He and his spouse loved that they could easily don/doff it on their own.”

Vicki Ralph, MPPA, OTR/L, CLT-LANA

Why choose ReadyWrap?

ReadyWrap features low-stretch straps that mimic short stretch bandaging — delivering effective compression therapy with high working and low resting pressure.

To support patient compliance, ReadyWrap offers color-coded VELCRO® straps for application and alternating strapping system to support one handed donning as needed.

With ReadyWrap, you can be certain your patients will have an effective and comfortable garment that easily adjusts to their limb and donning abilities.



Low-stretch material

Mimics short-stretch bandaging — the gold standard of treatment.



Color-coded straps

Provides easy step-by-step donning guidance.



50% Overlapping straps^s

Helps prevent gapping.¹



Blocked support spine

Helps maintain garment form to support all-day wear.¹




Foam Padding

Cushioning at elbow and wrist provides comfort and evenly distributes pressure.



Modular design

Garments work effectively on their own or in combination.



ReadyWrap Calf



ReadyWrap Arm

Indications

- Lymphoedema
- Chronic oedema
- Venous Disease

Available in:



A full holistic and vascular assessment should be carried out before treatment.

ReadyWrap

Adjustable compression system for upper limbs


Measuring Instructions

Circumference points:

- cA** widest width of the palm, just below the metacarpal phalangeal joints
- cC** narrowest point of the wrist
- cE** round the elbow joint and into the crease of the elbow
- cG** top of the arm, just below the axillary crease

Lengths:

- IA-C** the length between point cA and point cC
- IC-G** the length between point cC and point cG



ReadyWrap Arm

- Double-sewn liner allows for easy, one-handed donning and doffing
- Padded wrist and elbow to eliminate compression hot spots and provide additional comfort
- Bilateral design contours to both left and right arms

	Circumference (cm)			Length (cm) (IC-G)	REF	
	Wrist (cC)	Elbow (cE)	Axila (cG)		Beige	Black
Small	13.5-17	19-26	22-29	Short (up to 43cm)	104514	104515
				Average (up to 46cm)	104516	104517
				Tall (up to 49.5cm)	104518	104519
Medium	15.5-20	24-33	26-36	Short (up to 43cm)	104508	104509
				Average (up to 46cm)	104510	104511
				Tall (up to 49.5cm)	104512	104513
Large	17.5-23	29-40	31-43	Short (up to 43cm)	104502	104503
				Average (up to 46cm)	104504	104505
				Tall (up to 49.5cm)	104506	104507
X-Large	17.5-23	30-43	36-50	Short (up to 43cm)	104520	104521
				Average (up to 46cm)	104522	104523
				Tall (up to 49.5cm)	104524	104525



ReadyWrap Gauntlet

- Left and right specific design contours more naturally to the hand
- Built-in padding for palm and dorsum helps provide consistent compression over bony prominences
- Pocket design allows for foam pads to provide additional support
- Designed to overlap the ReadyWrap Arm garment to reduce gapping

	Circumference (cm)		Length (cm) (IA-C)		REF	
	Palm (cA)	Wrist (cC)			Beige	Black
Small	18-20	13.5-17.0	up to 12.5	Right	104498	104499
				Left	104490	104491
Medium	20-22	15.5-20.0	up to 13.5	Right	104496	104497
				Left	104488	104489
Large	22-24	17.5-23.0	up to 14.5	Right	104494	104495
				Left	104486	104487
X-Large	24-26	19.5-26.0	up to 15.5	Right	104500	104501
				Left	104492	104493



ReadyWrap

Adjustable compression system for lower limbs


Measuring Instructions

Circumference points:

- cG** 2 cm below Gluteal Fold
- cE¹** 5 cm above the Knee
- cD** 2-3 cm below Popliteal Fossa
- cC** widest part of Calf
- cB** 2 cm Above Malleolus
- cY** around Instep and Heel with foot dorsi-flexed
- cA¹** around the arch/dorsum of the foot
- cA** around base of Toes

Lengths:

- IA-Y** heel to base of Toes
- ID-B** the length between point cB and point cD
- IG-E¹** the length between point cE¹ (5 cm above the knee) and point cG



ReadyWrap Toes

- Unique design allows for adjustable toe girths
- Provides full foot coverage when paired with a ReadyWrap Foot unit
- Left and right specific design provides enhanced contouring

	Circumference (cm)		REF	
	Base of toes (cA)		Beige	Black
Small	18.0-20.4	Right	104416	104417
		Left	104406	104407
Medium	20.5-23.5	Right	104414	104415
		Left	104404	104405
Large	23.6-26.0	Right	104412	104413
		Left	104402	104403
X-Large	26.1-28.5	Right	104418	104419
		Left	104408	104409
XX-Large	28.6-31.0	Right	104420	104421
		Left	104410	104411



ReadyWrap Foot

- Left and right specific design and medial dart enhance contouring
- Easy to don, two strap design fits inside most shoes
- Overlaps with ReadyWrap Calf and ReadyWrap Toe units to help prevent gapping

	Circumference (cm)		Foot length (cm) (IA-Y)	REF		
	Anterior ankle/heel (cY)	Ankle (cB)		Beige	Black	
Small	28-33	18-23	up to 21	Right	104436	104437
				Left	104428	104429
Medium	30-38	20-28	up to 23	Right	104434	104435
				Left	104426	104427
Large	35-43	25-33	up to 26	Right	104432	104433
				Left	104424	104425
X-Large	40-48	30-38	up to 29	Right	104438	104439
				Left	104430	104431



ReadyWrap Foot SL

- Easy to fit in shoe with slim, tubular design
- Regular and Long lengths available to ensure coverage of the foot
- Left and right specific design to enhance contouring and improve fit
- Overlaps ReadyWrap Calf and ReadyWrap Toe units to help prevent gapping

	Circumference (cm)		Foot length (cm) (IA-Y)		REF	
	Foot (cA ¹)	Ankle (cB)			Beige	Black
Small	21-24	18-23	16-19cm (Regular)	Right	144425	144427
				Left	144424	144426
			19-22cm (Long)	Right	144429	144431
				Left	144428	144430
Medium	25-28	20-28	16-19cm (Regular)	Right	144433	144435
				Left	144432	144434
			19-22cm (Long)	Right	144437	144439
				Left	144436	144438
Large	29-32	25-33	16-19cm (Regular)	Right	144441	144443
				Left	144440	144442
			19-22cm (Long)	Right	144445	144447
				Left	144444	144446
X-Large	33-36	30-38	16-19cm (Regular)	Right	144449	144451
				Left	144448	144450
			19-22cm (Long)	Right	144453	144455
				Left	144452	144454



ReadyWrap Calf

- Can be worn with compression hosiery for additional support and protection
- Blocked spine reduces stretch and improves structural support
- Overlaps ReadyWrap Foot and ReadyWrap Knee garments to help prevent gapping

	Circumference (cm)		Length (cm) (IB-D)	REF	
	Ankle (cB)	Calf (cC)		Beige	Black
Small	18-23	28-38	up to 30	104450	104451
			up to 35	104452	104453
Medium	20-28	33-43	up to 30	104446	104447
			up to 35	104448	104449
Large	25-33	38-48	up to 30	104442	104443
			up to 35	104444	104445
X-Large	30-38	43-58	up to 30	104454	104455
			up to 35	104456	104457
XX-Large	35-43	53-68	up to 30	104784	104785
			up to 35	104786	104787



ReadyWrap Knee

- Easy, three-strap design provides a comfortable fit reducing pinching and high pressure points
- Overlaps ReadyWrap Calf and ReadyWrap Thigh to help prevent gapping

	Circumference (cm)		REF	
	Below Knee (cD)	Above Knee (cE ¹)	Beige	Black
Small	28-38	38-48	104464	104465
Medium	33-43	45-58	104462	104463
Large	38-48	52-68	104460	104461
X-Large	43-58	59-78	104466	104467
XX-Large	53-68	66-88	104788	104789



ReadyWrap Thigh

- Low stretch material provides low resting and high working pressures from superior patella border to the groin
- Proximal strap helps accommodate straight and conical thighs
- Overlaps with the ReadyWrap Knee to help prevent gapping

	Circumference (cm)		Length (cm) (IE ¹ - G)	REF	
	Above Knee (cE ¹)	Groin (cG)		Beige	Black
Small	38-48	48-63	up to 25	104478	104479
			up to 30	104480	104481
Medium	45-58	58-73	up to 25	104474	104475
			up to 30	104476	104477
Large	52-68	68-83	up to 25	104470	104471
			up to 30	104472	104473
X-Large	59-78	78-93	up to 25	104482	104483
			up to 30	104484	104485
XX-Large	66-88	88-103	up to 25	104790	104791
			up to 30	104792	104793



ReadyWrap Extender Strap

- Provides an additional 10cm to a standard ReadyWrap strap
- "Y" clamp design helps ensure a secure fit
- Compatible with ReadyWrap Calf, ReadyWrap Knee and ReadyWrap Thigh

	REF	
	Beige	Black
One Size	304398	304399



ReadyWrap Liners

Select the ReadyWrap Liner which is closest to the equivalent ReadyWrap garment

- Each ReadyWrap lower limb garment (excluding ReadyWrap Toes) is supplied with a ReadyWrap Liner
- Additional liners can be purchased in pairs
- Designed to be used under ReadyWrap garments to protect skin
- Provide no therapeutic compression

For ReadyWrap below knee garments

	REF
	Black
Small - Large	149032
X-Large - XX-Large	149033

For ReadyWrap Thigh garments

	REF
	Black
Small - Large	149034
X-Large - XX-Large	149035




Contact your ReadyWrap distributor for Measurement and Order Forms

Lower Limb

ReadyWrap™ Measurement and Order Form

How to measure the leg:
Ask the patient to stand, if possible. Mark and measure the leg as follows:



Write your measurement results in the box adjacent to measuring point

Measurement points required for each ReadyWrap garment	Circumference										Length		
	e1	e2	e3	e4	e5	e6	e7	e8	e9	e10	IA-C	IB-C	IC-G
ReadyWrap Ankle	x												
ReadyWrap Calf		x											
ReadyWrap Foot			x										
ReadyWrap Knee				x									
ReadyWrap Thigh					x								

Circumference points:
e3 2 cm below Gluteal Fold
e4 5 cm above the Knee
e5 2-3 cm below Popliteal Fossa
e6 widest part of Calf
e7 2 cm above Malleolus
e8 around instep and Heel with foot dorsiflexed
e9 around the arch/bottom of the foot
e10 around base of Toes

Lengths:
IA-Y Heel to base of Toes
IB-B The length between point e11 (2 cm above the malleolus) and point e12 (2 cm below popliteal fossa)
IC-G The length between point e11 (2 cm above the knee) and point e12 (2 cm below gluteal fold)


NB: e = circumference, l = length
All length measurements should be longitudinal and should follow body contour.
All measurements to be in centimetres (cm).

www.Lohmann-Rauscher.com.au

Upper Limb

ReadyWrap™ Measurement and Order Form

How to measure the arm:
Measure the arm in a supported and relaxed position with the elbow slightly flexed. The measurement points refer to the following anatomical landmarks:



Write your measurement results in the box adjacent to measuring point

Order Information

Order Date: _____ Patient Name: _____
Patient Phone: _____ Patient Email: _____
Clinician Name: _____ Institution/Organisation: _____
Clinician Phone: _____ Clinician Email: _____

Delivery

Attn: _____
Address: _____
City: _____ State: _____ Postcode: _____

Billing Information

Bill to: _____
Account #: _____
Card #: _____
Exp: / / COV: _____

Products

ReadyWrap Arm (per individual)	Colour	Quantity
<input type="checkbox"/> Black <input type="checkbox"/> Beige		

ReadyWrap Calf (per individual)	Style	Colour	Quantity
<input type="checkbox"/> Left <input type="checkbox"/> Right	<input type="checkbox"/> Black <input type="checkbox"/> Beige		

Lohmann Rauscher Pty Ltd, Suite 102, 21 Fernree Gully Road, Mount Waverley, Victoria 3149, Australia
Tel: +61 3 9578 8888 info@LohmannRauscher.com
www.Lohmann-Rauscher.com.au



To discover more about ReadyWrap and additional resources, visit www.SupportforLife.com

For all queries relating to ReadyWrap as well as SwellSpots, Caresia, TributeNight or Exo products please contact your local distributor:

Sources:

- Ehmann, S; Bock, K. (2017). Critical evaluation of bandage alternative characteristics: a case series. EWMA Conference, Amsterdam, Netherlands.
- Hampton, S. (2016). The difficulty and the solution of compression therapy in a healed venous leg ulcer. Br J Community Nursing, 21(9), p34-p38.
- Ehmann, S. (2016). Multinational, pilot audit of a VELCRO adjustable compression wrap system for venous and lymphatic conditions. JWC, 25(9), p513-p520.

* Lower extremity styles only.
§ ReadyWrap Calf, Knee and Thigh styles only.